




## EVALUATION OF DETECTABILITY, EXTENT, AND SEVERITY OF MYOCARDIAL INFARCTS IN PATIENTS

- **Receptor Quantification as a Function of Uptake Ratio(s)**
- **Partial Volume Effects**
- **Scatter and Attenuation – Correction Schemes**
- **Threshold for Changes in Uptake(s)**
- **Comparison of Different Acquisition Modes, e.g. 2D vs. 3D Pet(s)**

RSD's Heart & Thorax Phantom is designed to provide different reconstruction strategies, imaging protocol for patients, and testing and validation of image registration techniques. It also delivers a valid assessment of mammoscintigraphy techniques.



**Heart:** An accurately anatomic heart model is based on vacuum-formed shells. It was designed using high resolution, contrast-enhanced, ultrafast CT data from a normal patient, slightly modified to facilitate its use. The left and right chambers are connected at the atrium region to make a single compartment which can be filled and flushed independently using two ports labeled HC (heart chambers). The right ventricle is slightly modified to allow air to escape during filling. The myocardial wall (MW) has two ports, flushing and independent filling. The volume of the heart chambers is 284 ml, while the volume of the myocardial wall is 238 ml, without inserted defects. The standard model includes three defects with volumes of 8.9, 13.5, and 41.7 ml, respectively. Each of the defects can be filled separately. Defects of different dimensions can be ordered at no added cost. A disassembled heart is sent on request, so that dimensions of a special set can be established. Note that different defects cannot be retrofitted in the assembled heart.

**Basic Thorax:** The thorax is molded of polyurethane, modified for tissue equivalence, with a mass density of 1.10 g/cc. The narrow beam linear attenuation coefficient measured at 140 keV (Tc-99m) is 0.160 cm<sup>-1</sup>. The volume of the thoracic cavity, when all organs (heart, lungs, and liver) are inserted, is about 8,200 ml. It is filled from the top, with either an inert or a radioactive solution, to simulate the thoracic background. A valve is installed at the base for conveniently draining the phantom. The residue on the walls of the cavity can be easily flushed with the fittings provided at the top of the phantom. A second, smaller fitting is also provided as an air-bleed during filling.

**Lungs:** Perfusible lungs are molded in hollow, vacuum-formed shells, filled with Styrofoam beads. The final mass density of 0.40 g/cc is attained by adding an inert or radioactive solution through a filling port at the top of each lung shell. Extra sets of lungs can also be furnished for work continuity. The volumes of the left and right lung shells are 907 ml and 1,134 ml, respectively.

**Liver:** A liver with a volume of 980 ml is included to evaluate the effect of its uptake on quantitative myocardial imaging. It is a vacuum formed shell, mounted on acrylic tubes to minimize artifacts. The liver can be filled with an inert or radioactive solution.

**Fillable External Markers:** A set of fillable capsules is provided to serve as external markers. Capsules can be filled with a radioactive solution and attached to the external surface of the phantom. The phantom can then be imaged, using SPECT or PET modalities to compare image registration techniques.

**Thorax Overlay, Removable Breasts, and Breast Tumors:** The thoracic phantom without the overlay simulates an average male patient. The overlay, with or without breasts, simulates a large female or a still larger male patient, respectively. It is then possible to evaluate the effect of additional attenuation and scatter on quantitative myocardial imaging. The volume of each vacuum-formed breast is 972 ml. A tumor, filled with a radioactive solution can be inserted to evaluate the planar and tomographic imaging techniques used for mammoscintigraphy. A set of breast tumors (3, 6, 9, 12 and 15mm diameters) is included. They are supported by thin, threaded nylon rods which pass through ports and are sealed by O-rings. They can be bent by hand to reach any part of a breast.

## Model Numbers

Model No.	Product Description
RS-800T	Heart/Thorax Phantom
RS-801	Thoracic Cavity with the bottom plate
RS-802	Perfusable Lungs with Styrofoam beads (pair)
RS-804	Heart with two hollow defects in myocardial wall – standard size or to fit
RS-805	Liver Shell ONLY
RS-806	Chest Overlay
RS-807	Removable Breast with set of tumors

**Materials** See page 30 for more information.

RSD Soft Tissue	RSD Cortical Bone	RSD Trabecular Bone
•	•	•

## Specifications

Packing Size	Packing Weight
86W x 79D x 48H cm	32 kg
34W x 31D x 19H in	70 lb.

**Publication References:** 1) Gallivanone F, Carne I, Interlenghi M, et al. A Method for Manufacturing Oncological Phantoms for the Quantification of 18F-FDG PET and DW-MRI Studies. Contrast Media Mol Imaging. 2017;2017:3461684. Published 2017 Sep 7. DOI: <https://doi.org/10.1155/2017/3461684>. 2) Heikkinen JO, Kuikka JT, Rautio PJ. Interdepartmental audit with an anatomically realistic lung phantom. Journal of Nuclear Medicine Technology. 2006 Mar;34(1):34-42. PMID: <https://pubmed.ncbi.nlm.nih.gov/16517967>. 3) Doshi N, Basic M, Cherry S. Evaluation of the Detectability of Breast Cancer Lesions Using a Modified Anthropomorphic Phantom. Journal of Nuclear Medicine, 39, 1951-7. 1998. PMID: <https://pubmed.ncbi.nlm.nih.gov/9829588>.



## Applications


Myocardial perfusion SPECT image quality evaluation

Attenuation corrections

Evaluating patient doses

Training for SPECT/PET imaging

Quantitative phantom for SUV validation



## Modalities

SPECT/PET

# RSD

Radiology Support Devices, Inc.

# Thyroid Phantom



## Applications

I-131 therapy training

Thyroid uptake and scan quality assurance

SUV quantification and validation




## Modalities

SPECT/PET

## THYROID UPTAKE & IMAGE QUALITY ASSURANCE

The Thyroid Phantom is an integral part of the Fission-Product Phantom, but it can also be supplied separately (with or without a head). This phantom has the same contours as in the full Fission-Product Phantom and contains a neck with corresponding “Superhuman Skeleton” vertebrae, a hollow-shell thyroid with filling and flushing ports, and a front-cover of tissue-equivalent phantom material. This assembly has been designed for rapid removal of the cover plate and thyroid, a quick fill with an isotope solution and an equally rapid reassembly. This permits work with iodines of very short half-lives. A portion of the clavicles and sternum are included to further enhance the realism of the phantom.

**Publication Reference:** Kramer GH, Olender G, Vlahovich S, Hauck BM, Meyerhof DP. Comparison of the ANSI, RSD, KKH, and BRMD thyroid-neck phantoms for <sup>125</sup>I thyroid monitoring. Health Phys. 1996 Mar;70(3):425-9. DOI: <https://doi.org/10.1097/00004032-199603000-00014>. PMID: 8609037.



## Model Numbers

Model No.	Product Description
RS-542	Without Head
RS-545	With Head and Cervical Spine

**RSD Materials** See page 30 for more information.

RSD Soft Tissue	RSD Cortical Bone	RSD Trabecular Bone
•	•	•

## Specifications

Packing Size	Packing Weight
36W x 36D x 36H cm	5 kg
14W x 14D x 14H in	12 lb.

**Soft Tissues:** There are unlimited, small variations in density and absorption throughout the human body. Phantom soft tissue is closely controlled to have the average density of these tissues.

**Skeletons:** RSD skeletons are highly detailed polymer moldings which reproduce the shape, mass density and attenuation coefficients of cortical bone and spongiosa. RSD's proprietary moldings allow for continuous production, eliminate the restrictions of human skeleton bones (including limited availability, unethical collection of human bone specimen, variable size, and uncertain chemical composition), and avoid the loss of marrows in dried natural skeletons thereby making RSD skeletons superior to "real bone."

**Molds:** Molds for the RSD cortical bone and spongiosa were made from human skeletons consistent with the sizes of the soft tissue molds.

**ICRU 44:** RSD skeletons conform closely to the standards established by the International Commission on Radiation Units and Measurements ([ICRU Report No. 44](#)); mass density is reduced slightly to take into account a small decrease in calcium content for older patients.

## LINEAR ATTENUATION DATA

1. Monte Carlo simulation was used to calculate linear attenuation coefficients as a function of beam.
2. Monte Carlo results were validated with linear attenuation coefficients derived from Hounsfield Unit measurements at discrete energy levels.
3. RSD Phantom material linear attenuation data was compared to NIST data using ICRU Report 44 compositions of human tissues.
4. NIST data was interpolated when necessary.

MATERIALS	DENSITY (g/cc)
RSD Soft Tissue (Opaque)	1.08
RSD Soft Tissue (Transparent)	1.10
RSD Cortical Bone	1.18
RSD Trabecular Bone	1.17

RSD SOFT TISSUE					
Energy (MeV)	Mean (HU)	Calculated (M)	$\mu$ (ICRU 44)	% Difference	Ratio
00.08	60.30	0.1948	0.1932	0.0080	0.9921
00.10	52.88	0.1797	0.1795	0.0015	0.9985
00.12	57.10	0.1717	0.1709	0.0044	0.9956
00.14	52.95	0.1623	0.1624	0.0007	1.0007
00.20	--	0.1477	0.1439	0.0261	0.9746
00.30	--	0.1245	0.1246	0.0004	1.0004
00.60	--	0.0950	0.0941	0.0101	0.9900
00.80	--	0.0825	0.0826	0.0013	1.0013
01.00	--	0.0744	0.0743	0.0018	0.9982
02.00	--	0.0520	0.0519	0.0018	0.9982
03.00	--	0.0351	0.0357	0.0171	1.0174
06.00	--	0.0288	0.0291	0.0088	1.0088
08.00	--	0.0252	0.0255	0.0098	1.0099
10.00	--	0.0229	0.0232	0.0149	1.0151
15.00	--	0.0203	0.0203	0.0015	0.9985
20.00	--	0.0189	0.0189	0.0017	1.0017

RSD CORTICAL BONE					
Energy (MeV)	Mean (HU)	Calculated (M)	$\mu$ (ICRU 44)	% Difference	Ratio
00.08	1365	0.4345	0.4280	0.0151	0.9851
00.10	1048	0.3496	0.3562	0.0184	1.0188
00.12	0977	0.3211	0.3274	0.0191	1.0195
00.14	0902	0.2932	0.2986	0.0180	1.0184
00.20	--	0.2511	0.2513	0.0009	1.0009
00.30	--	0.2155	0.2137	0.0084	0.9916
00.60	--	0.1596	0.1598	0.0011	1.0011
00.80	--	0.1403	0.1402	0.0010	0.9990
01.00	--	0.1274	0.1261	0.0106	0.9895
02.00	--	0.0883	0.0885	0.0017	1.0017
03.00	--	0.0611	0.0625	0.0229	1.0235
06.00	--	0.0512	0.0525	0.0246	1.0253
08.00	--	0.0468	0.0474	0.0120	1.0121
10.00	--	0.0446	0.0444	0.0039	0.9962
15.00	--	0.0410	0.0409	0.0016	0.9984
20.00	--	0.0393	0.0397	0.0102	1.0103

RSD TRABECULAR BONE (SPONGIOSA)					
Energy (MeV)	Mean (HU)	Calculated (M)	$\mu$ (ICRU 44)	% Difference	Ratio
00.08	551	0.2849	--	--	--
00.10	515	0.2586	--	--	--
00.12	439	0.2337	--	--	--
00.14	318	0.1541	--	--	--