



RSD

Radiology Support Devices, Inc.

Fluoro/DR 1-Shot+



ADVANCED QA FOR FLUORO AND DR SYSTEMS IN JUST ONE SHOT

The Fluoro/DR 1-Shot+ is a highly versatile and reliable tool designed for routine quality assurance, ensuring consistent image quality and system performance. This phantom is ideal for testing spatial resolution, contrast resolution, image uniformity, and dose reproducibility in a variety of imaging applications. Whether you're using it for single-shot tests or ongoing system performance monitoring, the Fluoro/DR 1-Shot+ helps ensure your imaging systems are always operating to their full potential.

Product Features & Benefits

Compatible with Fluoroscopy & Digital Radiography: Ideal for quality assurance across multiple imaging systems

Comprehensive Evaluation: Supports tests for spatial resolution, contrast resolution, homogeneity, and image receptor dose

Improved Image Quality: Ensures consistent and high-quality diagnostic images through routine system evaluation

Robust Design: Built with durable materials to withstand repeated use while providing accurate, repeatable results

Easy Integration: Designed to seamlessly integrate into existing QA procedures, making it easy for imaging professionals to adopt and use effectively

Standards Compliance

IEC 61223-3-1: Evaluation and routine testing in medical imaging departments - Part 3-1: Acceptance tests - Imaging performance of X-ray equipment for radiographic and radioscopic systems

Specifications


- Dimensions: 180 mm (Diameter) x 10 mm
- 1 mm high purity Cu Attenuation Filter
- *Custom FOV sizes available upon request*

Accessories

- Comprehensive Product Manual
- Sturdy Carrying Case

Model Number

QA-FDR020: Fluoro/DR 1-Shot+



RSDphantoms.com

Call: 310-518-0527