

## MEETS CDRH TESTING REQUIREMENTS FOR SYSTEM SAFETY & PERFORMANCE

The Center for Devices and Radiological Health (CDRH) QA Phantom is a robust quality assurance tool designed to meet the demanding needs of radiography and fluoroscopy system testing. Engineered for precision and compliance with industry standards, this phantom supports comprehensive evaluation and maintenance of imaging systems.

### Product Features & Benefits

**Standard Compliance:** Aligns with CDRH testing requirements to ensure system safety and performance

**Versatile Design:** Suitable for a wide range of fluoroscopy and radiography systems

**Comprehensive QA Testing:** Facilitates the assessment of image quality, dose accuracy, spatial resolution, and contrast sensitivity

**Integrated Test Patterns:** Built-in patterns for evaluating field alignment, beam centering, and geometric distortion

**Optimized Imaging Performance:** Ensures systems deliver reliable, high-quality diagnostic images

**Cost-Effective Solution:** Detects system issues early to minimize downtime and repair expenses

### Model Number

**QA-FDR050:** CDRH QA Phantom

### Standards Compliance


**NEXT Protocol for 2003 Survey of Fluoroscopic X-Ray Systems** | IEC 61223-3-1: Evaluation and routine testing in medical imaging departments - Part 3-1: Acceptance tests - Imaging performance of X-ray equipment for radiographic and radioscopy systems | **AAPM Report No. 60:** Instrumentation Requirements of Diagnostic Radiological Physicists

### Specifications

- Dimensions: 180 mm x 180 mm x 10 mm
- Attenuation Filters
- Res-Mesh

### Accessories

- Comprehensive Product Manual
- Sturdy Carrying Case
- Fixed Probe Stand



**RSDphantoms.com**  
**Call: 310-518-0527**