



RSD

Radiology Support Devices, Inc.

Fluoro/DR 1-Shot



QA FOR FLUOROSCOPY & DIGITAL RADIOGRAPHY SYSTEMS IN JUST ONE SHOT

The Fluoro/DR 1-Shot is an essential tool for ensuring the quality and consistency of your fluoroscopy and tomography imaging systems. Designed for routine testing, it provides a comprehensive evaluation of image quality parameters, helping to identify performance variations and ensure optimal system functionality.

Product Features & Benefits

Comprehensive Evaluation: Allows testing of resolution, contrast, homogeneity, and dose accuracy in one device

Precision Indicators: Includes detailed markings for accurate measurement of system alignment

Streamlined QA: Simplifies routine quality assurance procedures, saving time and resources

Improved Diagnostic Accuracy: Ensures imaging systems deliver consistent, high-quality results, reducing diagnostic errors

Regulatory Compliance: Supports compliance with industry standards for radiography and fluoroscopy QA

Cost-Effective Solution: Minimizes downtime and maintenance costs by identifying performance issues early

Model Number

QA-FDR010: Fluoro/DR 1-Shot

Standards Compliance


IEC 61223-3-1: Evaluation and routine testing in medical imaging departments - Part 3-1: Acceptance tests - Imaging performance of X-ray equipment for radiographic and radioscopy systems | **DIN 6868-150:** Image quality assurance in diagnostic X-ray departments - Part 150: Acceptance test of medical radiographic and fluoroscopic X-ray equipment

Specifications

- Dimensions: 300 mm x 300 mm x 19 mm
- 0.05 mm tolerance
- *Custom FOV sizes available upon request*

Accessories

- Comprehensive Product Manual
- Sturdy Carrying Case



RSDphantoms.com

Call: 310-518-0527